



IITC animal restrainers are used with the noninvasive blood pressure systems for mice and rats they can also be used for tail vein injections or other studies where the animal needs to be restrained.

IITC carries seven different sizes; the weight range is from neonate rats and adult mice up to 1,000 gram rats. They are constructed of sturdy clear, 1/4" wall acrylic. Only a mild soap and water solution is to be used to clean the restrainers no ammonia based products should be used.

They are all supplied with an adjustable head gate and the end plate is removable by snaps. This allows the animal to "walk-in" as opposed to the traditional method of bootstrap holders where the animal is physically forced into the restrainer backwards. The backward forcing of animals may cause severe trauma.

Part #	Description	Size	Dimensions
84	Restrainer Adult mice and neonate rats	20 to 45 grams	1" ID (25 mm) - 4" Length
83	Restrainer Rodents	70 to 170 grams	1 1/2" ID (38 mm) - 6" Length
82	Restrainer Rodents	180 to 270 grams	2" ID (50 mm) -8" Length
81	Restrainer Rodents	250 to 400 grams	2 1/2" ID (63 mm) - 9" Length
80	Restrainer Rodents	400 to 600 grams	3"ID (75 mm) - 12" Length
79	Restrainer Rodents	600 to 800 grams	3 1/2" ID (88 mm) - 13" Length
78	Restrainer Rodents	750 to 1000 gram	4" ID (105 mm) - 13" Length

One year warranty